


rotator[®]

de-rotation straps

Rotator de-rotation straps are designed for pediatric patients with internal or external foot rotation that impairs gait. They are not recommended for patients with fixed rotational deformities.

Rotators are made with a latex-free material that combines comfort with the compressive elasticity needed to control rotation. Rotators can be worn next to the skin or overtop of clothing.

1 Fasten the waistband around the patient using the provided, double sided, hook pads and cut off any excess material.



2 Attach rotator straps to shoes by using the included D-rings. Attach as close to the toe as possible



3 Wrap the rotator straps around the legs to achieve the desired correction (fig. 4 & 5). Tension and number of wraps can be varied based on patient needs. (2 wraps above and below the knee are common).



A
For Internal Rotation Correction




B
For External Rotation Correction




LaunchPad
INNOVATING OUTCOMES

4

Secure the rotator straps to the waistband using the provided double sided hook pads and cut off any excess strap.



A
For Internal Rotation Correction



B
For External Rotation Correction



www.launchpad-europe.eu

Čeština :
Dansk :
Deutsch :
Eesti :
Español :
Français :
Íslenska :
Italiano :
Latviešu :
Lietuvių :
Magyar :
Norsk :
Nederlands :
Polski :
Português :
Română :
Slovenščina :
Slovenčina :
Suomi :
Svenska :
Türkçe :
Ελληνικά :
БЪЛГАРСКИ :
日本語 :



LaunchPad
INNOVATING OUTCOMES

www.launchpad-op.com



Manufactured by LaunchPad O&P
2545 Chicago Ave, Suite 412, Minneapolis, MN 55404 USA

Authorized Representative EU
Dirame Ortho, ASSESTEENWEG 27-29, B-1740 TERNAT, Belgium
T: +32 (0)2 582 82 50 / F: +32 (0)2 582 61 12

05-XINS01-02